

# MONITOR ARMS

## Monitor arms for diverse needs

With Merivaara monitor arms, equipment is easily adjustable to correct viewing and operating position and can be moved away if necessary. Soft adjustments and stable operation ensure that monitor arms remain in the desired position.

One or more displays can be installed on Merivaara monitor arms. All monitor arms are available as wall or ceiling mounted versions.





<b>Maximum load</b>	21 / 2x18 kg duo	30 kg	40 kg
<b>Weight ranges</b>	1-3,5 / 3,5-7 / 3-6 / 7-12 / 8-13 / 9-15 / 12-18 / 13,5-21 kg	7-11 / 8-12 / 12,5-18 / 17-26 / 23-30 kg	8-16 / 15-32 / 20-40 kg
<b>Lead-in cables</b>	Esim. 2 kpl max 9 mm	Esim. 1 kpl 9 mm + 2 kpl 3,5 mm + 1 kpl 6 mm ja 2 kpl 5 mm	Esim. 3 kpl max 12 mm + 2 kpl 6,5 mm
<b>Spring arm length</b>	910 mm	930 mm	988 mm
<b>Optional fixed arm length</b>	600 / 750 / 950 mm	700 / 850 / 1000 / 1150 mm	700 / 850 / 1000 / 1150 mm
<b>Additional information</b>	Yhdelle monitorille, saatavissa myös kahdella monitorivarrella	Yhdelle/kahdelle monitorille	Yhdelle/kahdelle monitorille
<b>Product Code</b>	Monitor carrier systems	Monitor carrier systems	Monitor carrier systems